

TDT Release 72 / 2.12

Released 1/3/2011

New in Release 72

Important!: As of release 72, drivers will no longer support the (Gigabit) PI5 interface.

- **RPvdsEx/Drivers:**
 - **Features:** Drivers have support for Optibit cards (PO5/PO5e) on all 64-bit Windows O/Ss. Support for the RZ6 Processor, RS4 Data Streamer, IZ2 Stimulator, RV2 Video Tracker.
 - **Components:** MCValList, MCLatch, MCSerSource, MCInt16ToFlt, MCInt8ToFlt, MCFloat2Int16, and MCFloat2Int8.
 - **Macros:** RZ6_AmpIn, RZ6_AudioIn, RZ6_AudioOut, RZ6_Control, Stream_Server_MC, IZ2_Control, Video_Access, RZ_Serial_Rec, RZ_Serial_Rec_noCmds, RZ2_Serial_Send, updated PZ2_Control and PZ3_Control for RZ_DSP_P support.
 - Support for PO5e/PO8e cards and LO5/LO8 interfaces.
 - New RVMap software included with TDT Drivers for configuring RV2.

New in Release 2.12

- **OpenEx:**
 - **New OpenBrowser Export Options:** OpenBrowser can now export directly to EDF and DDT format.
 - **OpenController Auto and Manual Save/Load States:** OpenController can now automatically load all of your settings as well as preserve them for future experiments.
- **OpenDeveloper:**
 - **64-bit Support:** OpenDeveloper and TTankMin support 64-bit TTankX.
- **OpenSorter:**
 - **Batch Statistics:** Statistics can now be batch calculated for a block with or without running a sorting algorithm. See the OpenSorter User Guide.
- **SpikePac:**
 - **Updates to Threshold Controls:** PCSort and BoxSort now offer improved threshold functionality and interfaces are more flexible.
 - **Control Set Save States:** Changes to PCSort and BoxSort control set settings can now be saved to OpenController save states.