


Fast Facts

RZ10 Lux Integrated Processor

This sheet provides basic information for the RZ10 and RZ10X Lux IO Processors and related devices. See the System 3 manual for more detailed information



RZ10(x) Hardware Connections

Using integrated LUX LEDs and Photosensors: Connect light output from LEDs directly to the appropriate fluorescent optical ports. Connect fluorescent response return directly to LUX PS1 photosensor. The FC connector has a small key that must be aligned to the cable. On LUX components this key is at the 10 o'clock position.

Using external devices:
Current output (M8 connector) – Connect directly to standalone LEDs

Voltage input/output (BNC) – Connect external photosensor to BNC input or LED driver to BNC output



Example Configuration: Fiber photometry 3 color, 1 subject/site

Other configurations may require a different setup

LUX LEDs & Photosensors

All LUX components for RZ10(x) driver and sensor slots. If needed, these can be changed out by unscrewing the module and removing the component with the RZ10 powered off.



Integrated LEDs



Integrated photosensor




Integrated Power Meter



4-pin M8 connector for standalone LED connection




BNC connector for external photosensor or LED driver connection



The RZ10(x) gizmo includes a LUX tab for automatically detecting and configuring all connected Lux components


Click 'Detect Hardware' and the RZ10(x) will autofill Row A and Row B with the appropriate component information




The RZ10(x) is typically used for Fiber Photometry applications.

Add the *Specialized* → *Fiber Photometry* gizmo to the RZ10(x) in the Processing Tree. Detected LEDs and sensors are enabled automatically

The RZ10(x) LUX configuration will automatically inform any gizmo targeting the LUX I/O



Driver signals are named by the LUX wavelength.



IMPORTANT! Assign the correct Lux I/O bank for each connected Fiber Photometry gizmo. The first connected gizmo assumes 'Upper Bank'; the second assumes 'Lower Bank'

Analog Input, ADC 4 and ADC 8



	Enable to...	Scaler	AutoID	ID	Api Acc
Adc.4	Single Chan	1	<input checked="" type="checkbox"/>	Adc4	
Adc.8	Off	1	<input checked="" type="checkbox"/>	Adc8	



Additional analog inputs are available via BNC ports ADC 4 (RZ10 and RZ10x) and ADC 8 (RZ10x only). These can be enabled in the RZ10(x) → ADC tab in Synapse


LUX Integrated PS1 Photosensor

Bandwidth (Hz)	DC - 700
Wavelength Range (nm)	320 - 1100
Gain	1x10 ¹⁰

LUX Integrated PM1 Power Meter

Bandwidth (Hz)	DC - 3000
Wavelength Range (nm)	320 - 1100
Gain	6.5 x10 ⁴

Digital I/O – Byte C, Bits 0 - 3



	Enable	Output	Invert	AutoID	ID	Epoc S
Port-A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PortA	Off
Port-B	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PortB	Off
Port-C.0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PortC0	Full


DB25 Digital I/O Connector Pinouts

Byte B			
Pin	Bit	Pin	Bit
10	1	22	0
11	3	23	2
12	5	24	4
13	7	25	6

Byte A			
Pin	Bit	Pin	Bit
6	1	18	0
7	3	19	2
8	5	20	4
9	7	21	6

Bit Addr Byte C			
Pin	Bit	Pin	Bit
14	1	1	0
15	3	2	2
16	5	3	4
17	7	4	6

Four bits of digital input/ output (I/O) can be accessed via BNC ports 0 – 3. All 24 bits of word-addressable or bit-addressable memory can be accessed via the 'Digital I/O' DB25 connector.



Spectral Response

