




Fast Facts

The Subject Interface Module (SIM)

This fast fact sheet provides basic information for the SIM. See the System 3 Manual and Synapse Manual for more



***Important** – Software configurations are intimately tied to the hardware setup of the device. SIM Sub Modules (IZV, PZD, PZA) will get **detected** in **Rig Editor** (menu → Edit Rig → Detect) of Synapse and gizmo options can be used accordingly.



Software Configuration

Use Sub Amps

Sub1
 Sub2
 Sub3
 Sub4
 Sub5
 Sub6
 Sub7
 Sub8

Enable

Channels: 16 Banks: A

Name ID: Auto Amp1

Options

Sampling Rate: System Rate 6 kHz

DC Coupled:

Reference Mode: Local

Filtering: 45% FS

External Ground:

Impedance Target: 1K

Set to base type: Single Unit

Use Sub Amps

Sub1
 Sub2
 Sub3
 Sub4
 Sub5
 Sub6
 Sub7
 Sub8

Enable

Boards: 1

Channels: 16 Banks: A

AC Coupled:

Name ID: Auto Amp1

Fast Access

Options

Sampling Rate: System Rate 6 kHz

Low Pass Filter: Auto

High Pass Filter: 0.1Hz




External Ground:

Impedance Target: 1K

See Synapse Manual for more details

PZA gizmo **PZD gizmo**

Stimulation Hardware/ Software Configurations



***Important** – Use Remote Arming to force user to arm manually from the touchscreen interface

See Synapse Manual for more details

SIM Stimulator Synapse Configurations. Use **Electrical Stim Driver** to send stim signals to the **IZV gizmo**. Set Normal, Repeated, Serial, and Parallel Mode in IZV gizmo. Serial and Parallel mode require special harnessing to achieve output.

SIM Stimulator/ Output Configurations for Normal, Repeated, Serial, and Parallel Mode. Serial and Parallel mode require special harnessing to achieve output.

Electrical Stimulation	
Stim Output Channels	16 per card
Stim Output Voices	4 per card
Stim Card Compliance	+/- 15 V* +/- 5 mA per voice, up to 3 kOhm load*
Stim Output Resolution	Voltage Mode: 100 uV Current Mode: 10 nA
DC Offset Current	Active channel: < 100 nA Open Channel < 1nA
Sample Rate	Up to 50 kHz
Analog Recording	
A/D	16 channels per card, hybrid A/D. Sample rate up to 50 kHz
Maximum Voltage In	+/- 500 mV
Frequency Response	DC Coupled: 0 Hz – 0.45*Fs (Fs = Sample Rate) AC Coupled: 0.4 Hz – 0.45*Fs
S/N (typical)	Single-Ended: 104 dB, Fs = 25 kHz, 300 Hz – 7000 Hz Differential: 116 dB, Fs = 750 Hz, 0.4 Hz – 300 Hz
Input Impedance	AC Coupled: 100 kOhm DC Coupled 20 MOhm
Digital Recording	
Input	Up to 128 Intan-based digital channels per bank
Maximum Voltage In	+/- 5 mV
Frequency Response	0.1 Hz – 10 kHz
Battery	32 Amp-Hour. 8-10 hours to charge to 95% capacity, 14 hours to full charge. <i>Battery life between charges dependent on # of active cards:</i> 2 active cards ~50 hr; 4 active cards ~35hr; 6 active cards ~27 hr; 8 active cards ~22 hrs
Charger	External 12 VDC, 2.5 A power supply, center negative

* Higher values can be obtained by combining voices and cards through harness and software configuration