

# Omnetics Analog Headstages

---

Hardware Reference



Updated 2025-02-27

© 2016-2025 Tucker-Davis Technologies, Inc. (TDT). All rights reserved.

Tucker-Davis Technologies  
11930 Research Circle  
Alachua, FL 32615 USA  
Phone: +1.386.462.9622  
Fax: +1.386.462.5365

### **Notices**

The information contained in this document is provided "as is," and is subject to being changed, without notice. TDT shall not be liable for errors or damages in connection with the furnishing, use, or performance of this document or of any information contained herein.

The latest versions of TDT documents are always online at <https://www.tdt.com/docs/>

# Table of Contents

---

**Omnetics Analog Headstages**

LP32CH-Z - 32 Channel Chronic Headstage	4
Headstage Voltage Range	4
Technical Specifications	5
RA16CH/LP16CH/LP16CH-ZNF - 16 Channel Chronic Headstage	6
Headstage Voltage Range	7
Technical Specifications	7

# Omnetics Analog Headstages

## LP32CH-Z - 32 Channel Chronic Headstage

The 32 channel chronic headstage is recommended for extracellular neurophysiology using silicon electrodes, metal microelectrodes or microwire arrays with input impedance from 20 kOhm to 5 MOhm. The headstage uses a female Omnetics connector to mate with chronic electrodes.

The LP32CH-Z uses an Omnetics 36 socket female dual row nano connector (0.025" or 0.64 mm) with 4 guide posts. The LP32CH-ZNF uses non-ferrous materials, including a non-ferrous female Omnetics connector with the same size and pinout as the LF32CH. It can be used for recording single-unit electrophysiology during fMRI neuroimaging.

Part Numbers:

LP32CH-Z - 32-Channel Chronic Low Profile Headstage for PZ PreAmps

LP32CH-ZNF - 32-Channel Non-Ferrous, Chronic Low Profile Headstage for PZ Preamps

### Warning

The headstage has sensitive electronics. Always ground yourself before handling.

## Headstage Voltage Range


When using a TDT preamplifier the voltage input range of the preamplifier (PZ5, Subject Interface, RA16PA) is typically lower than the headstage and must be considered the effective range of the system. Also keep in mind that the output range of the headstage varies depending on the power supply provided by the preamplifier. PZ5 and Subject Interface supply  $\pm 2.5$  V. PZ2 and RA16PA preamplifiers supply  $\pm 1.5$  V. Third party preamplifiers may vary. TDT recommends using preamplifiers which deliver  $\pm 2.5$  V or less. The table below lists the input voltage ranges for the LP32CH-Z headstage for either  $\pm 1.5$  V or  $\pm 2.5$  V power sources.

Power Source	Headstage Input Range
$\pm 1.5$ V	$\pm 1.48$ V
$\pm 2.5$ V	$\pm 2.49$ V

## Technical Specifications

Input referred noise	3 $\mu$ Vrms bandwidth 300-3000 Hz 6 $\mu$ Vrms bandwidth 30-8000 Hz
Headstage Gain	Unity (1x)
Frequency Response	DC - 25 kHz
Input Impedance	1e14 ohms
Connector	Omnetics 36 socket female dual row nano connector (0.025" or 0.64 mm) with 4 guide posts

## Pinout



The numbers on the pinout diagram above show the channel connections to the amplifier. By default, the headstage inputs are single ended, with Reference and Ground tied together by a jumper. To make the inputs referential, cut the jumper pictured below.

LP32CH:




### Important

When using multiple headstages, ensure that a single ground is used for all headstages. This will avoid unnecessary noise contamination in recordings. See the [Headstage Connection Guide](#) for more information.

## RA16CH/LP16CH/LP16CH-ZNF - 16 Channel Chronic Headstage

The 16 channel chronic headstages are recommended for extracellular neurophysiology using silicon electrodes, metal microelectrodes or microwire arrays with input impedance from 20 kOhm to 5 MOhm. They come in three configurations: standard profile, low profile, and non-ferrous low profile.

The RA16CH and LP16CH (low profile) use a female Omnetics 18 socket female dual row nano connector (0.025" or 0.64 mm) with 2 guide posts that is compatible with TDT microwire arrays, NeuroNexus chronic electrodes, and a wide variety of connectors. Users can also request a matching male Omnetics connector (OMCON\_ML\_HB) from TDT for custom built electrode arrays. The low profile LP16CH provides a smaller footprint than the RA16CH for better clearance in tight applications.



The LP16CH-ZNF uses non-ferrous materials, including a non-ferrous Omnetics connector with the same size and pinout as the RA16CH and LP16CH headstages. It can be used for recording single-unit electrophysiology during fMRI neuroimaging.

Part Numbers:

LP16CH - 16-Channel Chronic Low Profile Headstage for Medusa PreAmps

LP16CH-Z - 16-Channel Chronic Low Profile Headstage for PZ PreAmps

LP16CH-ZNF - 16-Channel Non-Ferrous, Chronic Low Profile Headstage for PZ Preamps

RA16CH - 16-Channel Chronic Headstage for Medusa PreAmps

RA16CH-Z - 16-Channel Chronic Headstage for PZ PreAmps

### Warning

The headstage has sensitive electronics. Always ground yourself before handling.

## Headstage Voltage Range


When using a TDT preamplifier the voltage input range of the preamplifier (PZ5, Subject Interface, RA16PA) is typically lower than the headstage and must be considered the effective range of the system. Also keep in mind that the output range of the headstage varies depending on the power supply provided by the preamplifier. PZ5 and Subject Interface supply  $\pm 2.5$  V. PZ2 and RA16PA preamplifiers supply  $\pm 1.5$  V. Third party preamplifiers may vary. TDT recommends using preamplifiers which deliver  $\pm 2.5$  V or less. The table below lists the input voltage ranges for the 16-channel chronic headstage for either  $\pm 1.5$  V or  $\pm 2.5$  V power sources.

Power Source	LP16CH Input Range	RA16CH Input Range
$\pm 1.5$ V	$\pm 1.48$ V	$\pm 0.9$ V
$\pm 2.5$ V	$\pm 2.49$ V	$\pm 1.9$ V

## Technical Specifications

Input referred noise	3 $\mu$ Vrms bandwidth 300-3000 Hz 6 $\mu$ Vrms bandwidth 30-8000 Hz
Headstage Gain	Unity (1x)
Frequency Response	DC - 25 kHz
Input Impedance	1e14 ohms
Connector	Omnetics 18 socket female dual row nano connector (0.025" or 0.64 mm) with 2 guide posts

## Pinout



The numbers on the pinout diagram above show the channel connections to the amplifier. By default, the headstage inputs are single ended, with Reference and Ground tied together by a jumper. To make the inputs referential, cut the jumper pictured below.

RA16CH:



LP16CH/LP16CH-ZNF:



#### Important

When using multiple headstages, ensure that a single ground is used for all headstages. This will avoid unnecessary noise contamination in recordings. See the [Headstage Connection Guide](#) for more information.