

P05/P05e/P05c/F05 Optibit Interface

Hardware Reference



Updated 2025-02-27

© 2016-2025 Tucker-Davis Technologies, Inc. (TDT). All rights reserved.

Tucker-Davis Technologies
11930 Research Circle
Alachua, FL 32615 USA
Phone: +1.386.462.9622
Fax: +1.386.462.5365

Notices

The information contained in this document is provided "as is," and is subject to being changed, without notice. TDT shall not be liable for errors or damages in connection with the furnishing, use, or performance of this document or of any information contained herein.

The latest versions of TDT documents are always online at <https://www.tdt.com/docs/>

Table of Contents

P05/P05e/P05c/F05 Optibit Interface

Optibit Overview	5
F05 Status LEDs	6
P05/P05e/P05c Technical Specifications	6

P05/P05e/P05c/F05 Optibit Interface



P05 Interface Card



P05e Interface Card



P05c Interface Card

Optibit Overview

The Optibit system (Optical Gigabit) is designed for users that require high-speed real-time control of System 3 devices or precise system-wide device synchronization. The Optibit interface consists of a PCI card (P05), PCIe card (P05e), or PCI Cluster card (P05c) that must be installed in the computer and one or more Optibit-to-zBus interface modules (F05) that mount in the rear slot of a zBus device chassis or is built into RZ Processors. When using the Optibit interface, all devices are automatically phase locked to a single clock.

Part Numbers:

P05 - Optical PCI Card for Optibit Interface


P05e - Optical PCI Express Card for Optibit Interface

P05c - Direct PCI Interface for Cluster Computing

F05 - Optibit to zBus Interface

Fiber Optic Connection

Devices are connected in a simple loop using provided high speed noise immune fiber optic cabling. See the [System 3 Install Guide](#) for installation instructions. For setting up the P05c card in Synapse, see [Cluster Processing with Synapse](#).



F05 Status LEDs


Four status LEDs on the face of the F05 indicate the connection status of the interface.

- | | |
|------------|---|
| Connected | The Connected LED is lit when the interface is powered on and the fiber optic cable labeled IN is connected properly. Although the Connected LED will light if only the IN cable is connected, both cables have to be connected properly for communication to take place. |
| Identified | The Identified LED lights when a software signal sent from the PC is recognized by the interface. This takes place when launching TDT software such as zBUSMon or running a Synapse experiment. |
| Activity | The Activity LED is lit when data is being sent to or from the TDT hardware. |
| Error | The Error LED lights when there is a connection or communication error. For example, this LED will light if the fiber optic cables are not connected properly. |

P05/P05e/P05c Technical Specifications

Interface transfer rates vary by transfer type and device. See [Interface Transfer Rates](#) for more information.

Below is a diagram of the compatible PCI and PCIe slots used with the P05 and P05e/P05c Optibit Interface cards.



P05

The P05 zBus to PC interface card must be installed in a standard size, compliant 3.3 V slot.

P05e / P05c

The P05e and P05c must be installed in a PCI Express slot. The P05e card uses a single lane (x1) but may be used in any PCIe slot size (x1, x2, x4, x8 or x16).




Important


Do not attempt to install in low-profile PCI slots. While low profile and standard PCI cards maintain the same electricals, protocols, PC signals, and software drivers as standard PCI expansion cards, the low profile bracket is not compatible with standard card size.

Fiber Optics

Standard cable length is 5 meters. Longer cables (up to 30 meters) are available on request.



P05e and P05c dimensions



P05 dimensions