

ZIF-Clip® Headstage Holders

Hardware Reference



Updated 2024-12-19

© 2016-2024 Tucker-Davis Technologies, Inc. (TDT). All rights reserved.

Tucker-Davis Technologies
11930 Research Circle
Alachua, FL 32615 USA
Phone: +1.386.462.9622
Fax: +1.386.462.5365

Notices

The information contained in this document is provided "as is," and is subject to being changed, without notice. TDT shall not be liable for errors or damages in connection with the furnishing, use, or performance of this document or of any information contained herein.

The latest versions of TDT documents are always online at <https://www.tdt.com/docs/>

Table of Contents

ZIF-Clip® Headstage Holders

Using the Holder with ZIF-Clip® Headstages	5
Using the ZCD-ROD32	7
Holder Dimensions	8
Z-ROD Dimensions (for analog and digital headstages)	8
ZCD-ROD Dimensions (for digital headstages)	9

ZIF-Clip® Headstage Holders



The ZIF-Clip® headstage holders securely hold your analog or digital ZIF-Clip® headstages during electrode insertion and can be used with most micromanipulators. The headstage holders, including the stabilizing rod, are approximately 4.5" in length. The stabilizing rod is 3" in length and has a 3/32" diameter. An aluminum lock pin ensures the ZIF-Clip® does not open during insertion.

Each holder is designed for use with the corresponding ZIF-Clip® or Digital ZIF-Clip® headstage.

Part Numbers:

Z-ROD32 - 16 or 32-channel analog ZIF-Clip® headstage holder (for ZC32)

Z-ROD64 - 64-channel analog ZIF-Clip® headstage holder (for ZC64 and ZD64)

Z-ROD96 - 96-channel analog ZIF-Clip® headstage holder (for ZC96 and ZD96)

Z-ROD128 - 128-channel analog ZIF-Clip® headstage holder (for ZC128)

ZCD-ROD32 - 32-channel digital ZIF-Clip® headstage holder (for ZCD32 and ZD32)

ZCD-ROD64 - 64-channel digital ZIF-Clip® headstage holder (for ZCD64)

ZCD-ROD96 - 96-channel digital ZIF-Clip® headstage holder (for ZCD96)


Using the Holder with ZIF-Clip® Headstages

Each holder is sized to fit a particular headstage and with the exception of the ZCD-ROD32 (see below), they all can be fitted to the headstage in the same way.


First, connect the probe or adapter to your ZIF-Clip® headstage **before** putting the headstage in the holder (the square guide provided to ensure the probe or adapter is connected with the correct polarity is hidden from view when the headstage is in the holder). See [ZIF-Clip® Headstage Adapters](#) for more information.

Next, gently slide the ZIF-Clip® headstage onto the holder until it is completely secure as shown in the images below.

Finally, secure the lock pin to the headstage holder.



Gently slide the headstage onto the holder (with probe or adapter already connected).



Position the headstage holder between the cables of the ZIF-Clip® headstage. The headstage should be completely secured in the holder.




To remove, grip the top and bottom of the headstage and gently slide the holder away.



The U-shaped lock pin secures the connection and prevents the ZIF-Clip® from opening and releasing the probe.

Using the ZCD-ROD32

The ZCD-ROD32 has a unique design that requires a different insertion procedure.




To use the headstage holder:

1. Set the ZCD32 headstage inside the base (or U) of the holder and slide it forward until it is stopped by the interior flange (Image 1-4).
2. After the clip is in place, insert the probe (Image 5-6) and then slide the provided lock pin over the ZCD32 (Image 7-9).

The lock pin prevents the clip from opening and releasing the probe, and also from sliding backward during insertion.

The lock pin has small ridges that should be aligned with the grooves on the face of the clip. If you have trouble connecting the lock pin, make sure that the clip has been pushed in completely and that the ridges and grooves are properly aligned (Image 7).

Holder Dimensions



Z-ROD Dimensions (for analog and digital headstages)

Form Factor	16/32-channel ZC16 / ZC32	64-channel ZC64 / ZD64	96-channel ZC96 / ZD96	128-channel ZC128
Height	4.10 mm (9.62 mm with lock pin)			
Inner Width	9 mm	14 mm	17.50 mm	24 mm
Outer Width	13 mm (16 mm with lock pin)	18 mm	21.50 mm	28 mm
Holder Length	25 mm	28 mm	28 mm	28 mm
Rod Length	76.2 mm			
Rod Diameter	2.29 mm			

ZCD-ROD Dimensions (for digital headstages)

Form Factor	32-channel ZCD32 / ZD32	64-channel ZCD64	96-channel ZCD96
Height	5.5 mm (11.50 mm with lock pin)	4.10 mm (15.30 mm with lock pin)	4.10 mm
Inner Width	11.10 mm	14.39 mm	17.50 mm
Outer Width	18.50 mm (18.50 mm with lock pin)	17.79 mm (21.5 mm with lock pin)	21.50 mm
Holder Length	31.50 mm	25.36 mm	25.36 mm
Rod Length	76.2 mm		
Rod Diameter	2.29 mm		