

# UZ3 USB 3.0 Interface for Optibit



UZ3 Interface

## UZ3 Overview

The UZ3 is a fiber optic interface to a high speed USB 3.0 port on your laptop or PC. It can connect to one or more RZ Processors. RX and RP devices are not supported by this interface.

A USB Type C-C cable and C-A adapter are included.

Part Numbers:

UZ3—USB 3.0 Interface

## UZ3 Technical Specifications

Interface transfer rates vary by transfer type and device. See “Cycle Usage and Large Transfers” on page 20-4, for more information.

The UZ3 can handle roughly 50-75% of the data rate of the PO5/PO5e card and should not be used for high channel count applications. The actual maximum rate is dependent on your PC specs.

The status LED is on by default when connected and flashes off when there is transfer activity.

Note: In UZ3 serial number 1021 and below, the LED is off by default and flashes on when there is activity.

Do not use with a USB 2.0 port. Use zBusMon software to measure the transfer rate to see if you are getting high speed data transfer.

The UZ3 dimensions are 31mm x 21mm x 106mm.

