


# LXFX-KIT-2C

Two Color Fiber Photometry Cable Kit



**Drive Cables:** 200um x 2m (0.5NA)  
LEDs => Optical Filter



**Response Cable:** 600um x 2m (0.5NA)  
Optical Filter => Photosensor




**Subject Cable:** 400um x 2m LAF (0.5NA)  
Optical Filter => Subject/Power Meter



**Cable Coupler:**



**Attenuating Coupler:**  
~90% with 200um cable




# LXFX-KIT-3C


Three Color Fiber Photometry Cable Kit



**Drive Cables:** 200um x 2m (0.5NA)  
LEDs => Optical Filter




**Response Cables:** 600um x 2m (0.5NA)  
Optical Filter => Photosensors



**Subject Cable:** 400um x 2m LAF (0.5NA)  
Optical Filter => Subject/Power Meter



**Cable Couplers:**



**Attenuating Couplers:**  
~90% with 200um cables

