

Self-Administration

Description	TDT iCon Module				Compatible Equipment		
	iHx	iR5	iMx	iL24	Lafayette	Med Associates	Harvard Apparatus / Coulbourn
Lights (Task)	X			X	80226xxx	ENV-321M	LE100567
Lever (Response)	X			X	80110	ENV-310W	LE100265
Pellet Dispenser (Reward)	X				80209-xxx	ENV-203-20	LE100250
Syringe Pump (Reward)	X					PHM-100VS-2	

Forced Choice

Description	TDT iCon Module				Compatible Equipment		
	iHx	iR5	iMx	iL24	Lafayette	Med Associates	Harvard Apparatus / Coulbourn
Lights (Task)	X			X	80226xxx	ENV-321M	LE100567
Lever (Response)	X			X	80110	ENV-310W	LE100265
Pellet Dispenser (Reward)	X				80209-xxx	ENV-203-20	LE100250

Acoustic Startle

Description	TDT iCon Module				Compatible Equipment		
	iH10	iR5	iMx	iL24	Lafayette	Med Associates	Harvard Apparatus / Coulbourn
Speaker (Startle Stimulus)			X		80135Mxxx	ENV-223, ENV224	
Force plate transducer and amplifier (Response)			X		HSM1000	PHM-250	LE 111

Fear Conditioning

Description	TDT iCon Module				Compatible Equipment
	iHx	iR5	iMx	iS9	Med Associates
Speaker (Neutral Stimulus)			X		ENV-223, ENV224
Electrical shocker (Aversive Stimulus)	X			X	NV-410C, ENV-413C
Force plate transducer and amplifier (Response)			X		PHM-250